

Verified by: Jyothi Vishumangalam Quality Release date: May 25, 2017

J. Vishminangalam.

Chai Biotechnologies Inc. Santa Clara, CA 95054 Phone: +1 (650) 779-5577

Product Name: FAM + HEX Calibrator Set

Catalog Num: S0140 Lot Num: 4079923

Analysis Date: May 25, 2017 Expiration Date: May 19, 2019

Storage: Store at 4 to 20 °C before reconstitution; 4 °C after

reconstitution; protect from light

Functional Test Results Fluorescence Spectrum of FAM Calibrator The freeze-dried FAM calibrator was reconstituted in the kit buffer as per instructions in the manual and a fluorescence emission spectrum was measured at 1 nm intervals, at a scan speed of 20 nm/min using 470 nm **PASS** excitation wavelength. The peak emission wavelength is specified to be $522 \pm 5 \text{ nm}$. Fluorescence Spectrum of HEX Calibrator The freeze-dried HEX calibrator is reconstituted in the kit buffer as per instructions in the manual and a fluorescence emission spectrum is measured at 1 nm intervals, at a scan speed of 20 nm/min. Using 470 nm **PASS** excitation wavelength. The peak emission wavelength is specified to be $557 \pm 5 \text{ nm}$.