

Verified by: Jyothi Vishumangalam

J. Vishnunangalan.

Quality Release date: Oct 27, 2017

Chai Biotechnologies Inc.

Santa Clara, CA

Phone: +1 (650) 779-5577

Product Name: Fluorescein Calibrator

Catalog Num: S0110

Lot Num: 8501669

Analysis Date: Oct 27, 2017

Expiration Date: Oct 27, 2019

Storage: Store at -20 °C

Functional Test Results

## Fluorescence Spectrum of Fluorescein Calibrator

The Fluorescence Emission spectrum of the Fluorescein Calibrator was measured at 1 nm intervals, at a scan speed of 20 nm/min using 490 nm excitation wavelength. The peak emission wavelength is specified to be  $521 \pm 5$  nm.

**PASS**